

## Weigh Dishes with Pour Spout



Weigh Boat Kit Anti-Static Polystyrene 30/pkg (10 x each of SM,MD,LG)	YBW01K
Small Anti-Static Polystyrene Weigh Boat- 50/pkg (.8125" x 2.25" x .5" )	YBW01S
Medium Anti-Static Polystyrene Weigh Boat - 50/pkg (3.25" x 5.125" x 1.25")	YBW01M
Large Anti-Static Polystyrene Weigh Boat - 50/pkg (4.25" x 7" x 1.125")	YBW01L

## Printers

Dot matrix impact printer provides long-lasting records. Includes interface cable and 1 roll of paper.

Dot-Matrix Printer for Pinnacle Series Balances	YDP01P1
Dot-Matrix Printer for Summit and XP Series Balances	YDP01XP
Dot-Matrix Printer for Timberline Series Balances	YDP05
Dot-Matrix Printer for MAXX Series (includes interface)	YDP02MXX
Dot-Matrix Printer for Mesa Series Balances (cable not included)	YDP02
GLP Printer with Time/Date Stamp	301380.1
Paper for Dot Matrix Printer (1 roll)	901044.1
Ribbon for Dot Matrix Printer (1 roll)	901045.1



## Draft Shields


Glass cabinet draft shield for use with Pinnacle and Summit Series mg balances. Housing with glass doors on top and sides ensures stable readings.



Pinnacle Series Balances Standard Draft Shield (mg models) 6.1" x 9.4" x 8.0" (156 x 238 x 203mm)	YDS03PI
Pinnacle Series Balances Mini-Draft Shield 6.1" x 6.5" x 8.0" (156 x 238 x 229mm)	YDS01PI
Pinnacle Series Balances (older mg models) 6.1" x 9.4" x 9.0" (156 x 238 x 229mm)	602643.1
Summit Series Balances (mg models) 6.5" x 6.6" x 9.5" (165 x 168 x 241mm)	YDS01S

## Draft Rings

Acrylic ring with stainless steel cover ensures stable readings. Fits round pan models only.



Pinnacle Series Balances 4.5" dia. x 1" height (114 x 25mm)	YDS02PI
Pinnacle Series Balances (older models) 6.0" dia. x 3.0" (152 dia. x 60 mm)	902191.1
Summit Series Balances 6.4" x 2.4" (161.5 x 60 mm)	YDS02S
MAXX Series Balances (mg units only) 4.4" dia. x 1" height (122 x 25mm)	69V00003

## Under Counter-Mounted Anti-Theft Lockdown Kits

Secures Denver balances to benchtop. Includes two hardened-steel cylinders, two keys to lock the balance in place, and template for drilling.

Pinnacle and XP Series Balances	36800110.1
---------------------------------	------------

## Anti-Theft Cable and Lock (not shown)

Secures balances to benchtop. Thread cable through loop on balance and around a sturdy object, and lock.

Pinnacle and XP Series Balances	400171.1
Timberline Series Balances	U1-LD2





### Weigh-Below Kit

Used for suspending samples below the balance to determine specific gravity. Includes 75mm threaded rod with hook. Connects directly to weigh below port.

Pinnacle Series Balances	YHW01
--------------------------	-------

### Density Determination Kit

Measures the density of solid, liquid or powder samples. Attaches directly to the balance for the most precise results.

Pinnacle Analytical Balances	301377.1
Summit Analytical Balances	YDK01S



### In-Use Covers

Durable, high-impact plastic cover to protect your balance during use.

Pinnacle Series Balances Analytical or Round Pan	602619.1
Pinnacle Series Balances with Square Pan	602620.1
Summit Series Analytical Balances	YDC04S
Summit Series mg, Round Pan	YDC03S
Summit Series (402 and 602 models)	YDC01S
Summit Series mg, Square Pan	YDC02S
Timberline Series Toploading Balances w/Square Pan	6960BL02
Timberline Series Toploading Balances w/Round Pan	6960BL03

### Serial Interface Cables

All cables are 6' (1.8 meters) in length and interface to computers with 9-Pin connectors.

Pinnacle Series Balances	YCC99-0006
Summit and XP Series Balances	YCC01-SM1
Timberline Series Balances	U1-LAT-CAB
MAXX Series RS232 (port w/cable)	YADAP-RS
MAXX Series USB (port w/cable)	YADAP-USB



### Labtronics Collect Software

Connect any Denver Instrument balance, with RS232, to your PC for direct input of real-time data into spreadsheet format (Interface cable not included).

Collect 6.1	902242.1
Collect SL	YSW04-0002

### Anti-Static Brush

Eliminates dust particles and neutralizes statically charged environments. Ionize the surrounding air and avoid unwanted foreign particles from introducing static interferences.

Anti-Static Brush for all Balances	69US0094
------------------------------------	----------



### Anti-Vibration Granite Isolators

Durable with a small footprint. Includes a granite top, silicone gel isolators, and ABS base construction. Slides under the balance to level any surface and reduce vibration interference.

13" x 15" Anti-Vibration Platform	U1-2120315
16" x 21" Anti-Vibration Platform	U1-24201621